





# Rule out canine KCS with the Schirmer Tear Test



The Schirmer Tear Test is comprised of a paper strip that is folded and hooked over the lower eyelid. Once in place, tears begin to “wick” up the filter paper. The distance, in millimeters, that the tears wick up the test strip in one minute, represents the millimeters of tear production.

## Early diagnosis improves treatment outcome

STT Result	Interpretation	Follow up*
$\geq 15$ mm/minute 	Normal	<ul style="list-style-type: none"> <li>No clinical signs - Retest at next exam</li> <li>Clinical signs - Further diagnostics</li> </ul>
$> 10-14$ mm/minute 	Early or sub-clinical KCS	<ul style="list-style-type: none"> <li>No clinical signs - Retest in 4 weeks</li> <li>Clinical signs - Start Treatment</li> <li>Exophthalmic breeds - Start Treatment</li> <li>Decreasing tear production has been documented - Start Treatment</li> </ul>
$6-10$ mm/minute 	Mild to moderate KCS	Start treatment
$\leq 5$ mm/minute 	Severe KCS	Start treatment

Although Keratoconjunctivitis sicca (KCS) is bilateral, each eye may be at a different stage in the disease process at the time of examination. Therefore, if one eye tests abnormal, the other eye should also be treated or monitored closely.

\* The above graphic illustration is a general guide to Schirmer Tear Test interpretations and suggested follow-ups. It is recommended veterinarians make their own assessment for each individual case.

## Progression of Clinical Signs

### Healthy Eye



#### Healthy Eye

- Normal tear production  
STT >15 mm

### Early KCS



#### Early KCS

- Slight reduction in tear production, STT <15 mm
- Possible minor irritation
- Vision not, or only slightly, affected
- Signs are not always noticeable
- Increased susceptibility to infection
- Treatment can resolve signs and treat the underlying cause

### Intermediate KCS



#### Intermediate KCS

- Insufficient tear production, STT <10 mm
- Irritation of the eye
- Loss of vision possible
- Frequent recurrent infections
- At risk of ulceration, often slow to heal
- Treatment can stabilise the disease and prevent recurrent infections

### Advanced KCS



#### Advanced KCS

- Severe reduction in tear production, STT <5 mm
- Severe irritation of the eye and formation of crusting
- Reduction of vision
- Corneal oedema, neovascularisation and pigmentation
- High risk of ulceration, slow to heal, can progress to perforation
- Treatment partially restores sight and prevents further loss of vision

## Schirmer Tear Test

There are 10 envelopes per box, each envelope has 5 sets of 2 strips, each set of 2 in a sterile packet. Each strip is calibrated in mm for easy measuring, has right and left eye markers and a blue dye that starts wicking at 5mm for easy tracking of the measurement.



\* NADA 141-052, Approved by FDA

Copyright © 2011 Intervet Inc., a subsidiary of Merck & Co., Inc. All rights reserved.  
Intervet Inc. d/b/a Merck Animal Health, Summit, NJ 07901 MAH-OPT-09